

Fetal Monitor



Cost-Effective Single Fetus Monitor

FC700

- FHR(Fetal Heart Rate), UC(Uterine Contraction) and FM(Fetal Movement)
- Automatic Fetal movement detection
- LED display of FHR and UC
- Alarm on abnormal FHR
- Easy-to-read A4 Size chart paper
- Use inexpensive fax paper

- ✓ High Sensitivity(1 MHz)
- ✓ S/W upgrade through RS-232C
- ✓ Auto CTG Analysis
- ✓ Wireless central monitoring up to 16 beds (Option)



Bionet
www.ebionet.com

SPECIFICATIONS

■ General Specifications																	
Dimension	296(W) × 305.5(H) × 92.5(D)mm, Approx.2.9kg																
Display	7 -Segment LED 2 channels(FHR, UC)																
Recorder	Method : Thermal print Print speed : 1, 2, 3 cm/min Paper feeding function																
■ Power Specifications																	
Power Adaptor	Input : 100 - 240VAC, 50 - 60Hz, 1.2A Output : 18VDC, 2.5A																
External Link	RS-232C : Program Download, Central System Connection(option)																
■ Environmental Specifications																	
Temperature Range	Operating temperature : 10 ~ 40°C Storage temperature : -10 ~ 60°C																
■ Performance Specifications																	
Fetal Heart Rate	Input signal : Ultrasound pulsed doppler Ultrasound frequency : 1.0Mhz Ultrasound power : <10mW/cm ² FHR detection method : Auto correlation FHR range : 50 - 240 bpm FHR accuracy · 120 - 160 bpm : ±1 bpm · Except 120 - 160 : ±2 bpm																
Uterine Contraction	Input source : External transducer with strain gauge Reference(zero) control : One touch switch Measurement range : 0 - 99 units																
Fetal Movement Measurement	Detection source : Ultrasound pulsed doppler																
Auto CTG Analysis	Average Baseline FHR Number of UC Number of Acceleration Number of Deceleration · Late Deceleration · Early Deceleration · Variable Deceleration Tachycardia(Moderate, Severe) Bradycardia(Moderate, Severe) *CTG Analysis results is printed out every 10 minutes. (Intermediate Report)																
■ Accessory																	
Standard	<table border="0"> <tr> <td>Ultrasound doppler probe</td> <td>1ea</td> </tr> <tr> <td>UC probe</td> <td>1ea</td> </tr> <tr> <td>Event marker</td> <td>1ea</td> </tr> <tr> <td>Print paper</td> <td>2rolls</td> </tr> <tr> <td>Power adaptor & cord</td> <td>each 1ea</td> </tr> <tr> <td>Ultrasound gel</td> <td>1ea</td> </tr> <tr> <td>Probe belt</td> <td>2ea</td> </tr> <tr> <td>Operation manual</td> <td>1ea</td> </tr> </table>	Ultrasound doppler probe	1ea	UC probe	1ea	Event marker	1ea	Print paper	2rolls	Power adaptor & cord	each 1ea	Ultrasound gel	1ea	Probe belt	2ea	Operation manual	1ea
Ultrasound doppler probe	1ea																
UC probe	1ea																
Event marker	1ea																
Print paper	2rolls																
Power adaptor & cord	each 1ea																
Ultrasound gel	1ea																
Probe belt	2ea																
Operation manual	1ea																



*Specifications subject to change without prior notice.

*You may have distortion or signal noise when you use nonstandard or other brand's accessories.

We strongly recommend you use only the authorized accessories which we supply.

شرکت ارکان آرا تجارت البرز

Bionet

Iran : Arkan Ara Tejarat Alborz Co., Ltd.

Tel : 021-81060
E-mail: info@arkanara.com
www.arkanara.com

South Korea : Bionet co., Ltd.

Tel : +82-70-7585-6418(Int'l) / +82-2-6300-6483(Domestic)
E-mail : sales@ebionet.com(Int'l) / yunms@ebionet.com(Domestic)
www.ebionet.com

North America(U.S.A) : Bionet America, Inc.

Tel : +1-714-734-1760
E-mail : sales@BionetUS.com
www.BionetUS.com

Latin America : Bionet Latin America, Inc.

Tel : +1-949-800-9585
E-mail : jayshim@ebionet.com
www.ebionet.com